

## DAEYOMEDI CO., LTD.

DAEYOMEDI Co., Ltd. is a company developing professional medical device and we are concentrating our efforts to diagnosis cardiac health and oriental medicine pulse diagnosis by using 3-dimensional blood pulse measuring and analysis technology. Vision of DAEYOMEDI Co., Ltd. is Worldwide leader of medical service company with respect to human life through technology and confidence.



[DMP-3000]



3-D Blood Pressure Pulse Analyzer

3-D Blood Pressure Pulse Analyzer 3-D Blood Pressure Pulse Analyzer

Specification	
Product name	DMP-Life
Medical Class	Class II-a
MDD classification	Pulse waveform analyzer
Power	AC110~220V
Range of Pressure	0~400g.f
Range of beat	40~200bpm
Accuracy of beat	5%
Accuracy of pressure	5%
Communication	USB to PC
Principle of measurement	Tonometry
Size(mm)	420 X 250 X100
Weight	4Kg







#### DAEYOMEDI CO., LTD.

-- 0 -- 0

Gyeonggi Technopark Technology Development Center RM601. 705 Haean-ro(st), Sangnok-gu, Ansan-si, 426-901, Korea http://www.daeyomedi.com dadaeyo@daeyomedi.com









## **DMP-Life is**

#### Swiss made accurate robotic cuff system

- Precision mini robot module can take your pulse more stably

#### Very Efficient Monitoring system

- DAEYOMEDI's patented technology make you more quick, more easy and more clear. - DMP-Life will save your money, time and human resource. You need just 1 system for CVD screening.

## **DMP-Life for**

## Hemodynamic Monitoring & Screening - DMP-Life can measure not only blood pressure but also vessel's condition, systemic circulation resistance

- and pulse condition.
- You can use results from DMP-Life to diagnosis patient's general health state

#### Health Promotion Center

- Sports Center, Spa and Healing center also can use DMP-Life - Health promotion center can use results from DMP-Life to predict CVD risk rate index(in case of use of
- Framingham score) or make regimen

# **DMP-Life has 4 function in 1 system**

## **Pulse Diagnosis**

## STRENGTH

**Powerful** : Increase of pulse pressure or increase of stroke volume or Inflammation. Increase of Heart load

**Powerless** : Decrease of pulse pressure or decrease of stroke volume. Lowering of heart pumping function

## DEPTH

Shallow : Fever or inflammation Acute course

**Deep** : Sleepless or chronic disease. Lowering of recoverative powers

## HEART RATE

Fast : Inflammation, fever or pyrexia. Increase of heart work

**Slow** : Fatigue, Cold symptom, aging or heart disease

## SHAPE

**Slippery** : Increase of stroke volume or blood plasma

**Choppy** : Anemia, Decrease of blood plasma and pulse pressure



#### **ECRI** Estimated Circulation Resistance from Flow/Pressure

Item	Meaning	Normal Range
ESV	Estimated Stroke Volume	60~100 ml/beat
SVI	Stroke Volume Index	33~47ml/beat/m^2
ECO	Estimated Cardio Output	4~8 L/min
ECRI	ECR Index	1,600~2,560 dyne*sec*m^2*cm^5

- Function of the circulation is transport blood to deliver oxygen, nutrients and chemicals to the cells of the body to ensure their survival and proper function and to remove the cellular wastes
- Resistance is inversely proportional to radius, resis tance is also affected by the viscosity of the blood flowing

DMP-Life realizes tonometry blood pressure technology(most precision method to measure blood pressure) by using array pressure sensor and precision robot system

- For maintaining health state every normal healthy person need to care there cardiovascular condition



## **Blood vessel**



Ref. Journal of Korean Medicine 2009

- If blood vessel looses elasticity reflected wave returns guickly and wave peak is more higher
- Stress and mental or social pressure, sleepless, smoking and aging can be a main cause of vessel stiffness.

## **Blood pressure**



[Tonometry method to measure Blood Pressure]